



**SuperVap<sup>®</sup> PFC**  
Concentration System for 15 and 50ml Centrifuge tubes  
Designed for PFAS Analysis

# The SuperVap® PFC Family of Concentrators

Designed specifically for the concentration and evaporation of samples for PFAS analysis with no components that will contribute to PFAS background. The SuperVap® PFC 24 Concentrator is a dry, waterless system capable of automatically concentrating or evaporating up to 24 samples in 15ml centrifuge tubes. The SuperVap® PFC 12 Concentrator can automatically concentrate or evaporate up to 12 samples in 50 ml centrifuge tubes. It can preheat as well as ramp up to the final temperature. It automatically starts the Nitrogen blowdown and shuts off Nitrogen and heat when the final programmed time is achieved. Samples can be concentrated directly to the centrifuge for unattended transfer, eliminating contamination and errors during manual transfer.

## **Dry, Waterless System:**

The SuperVap® PFC Concentrator operates without water, ensuring a dry environment for sample concentration.

## **Programmable:**

Users can program the device to perform specific tasks and follow a set sequence of steps. This flexibility allows for customization based on the specific requirements of the samples being processed. All Steps are programmed by time.

## **Preheating and Temperature Ramp-Up:**

The concentrator can preheat and gradually ramp up to the final temperature needed for the concentration process. This feature may be helpful for sensitive samples that require controlled temperature conditions.

## **Automated Nitrogen Blow Down:**

The device automatically initiates the Nitrogen blow-down process. Nitrogen is commonly used to evaporate solvents from samples during concentration.

## **Automatic Shut-Off:**

The concentrator is programmed by time to automatically shut off the Nitrogen and heat when the final programmed time is reached. This ensures the concentration process is precisely controlled and does not exceed the specified duration.

## **Unattended Transfer to Centrifuge:**

After concentration, the samples can be transferred directly to a centrifuge without manual intervention. This feature reduces the risk of contamination and minimizes errors during manual transfer processes.

## **Error Reduction:**

The automation and programmability of the SuperVap® PFC Concentrator contribute to the reduction of contamination and errors associated with manual sample transfer. This is particularly important in maintaining the integrity of experimental results.

## **Contamination Reduction:**

The SuperVap® PFC Concentrator reduces contamination. It is constructed with parts that do not contribute to the PFAS background. The materials are Peek®, Delrin®, LDPE. This is particularly important in maintaining the integrity of experimental results. Hepa and Carbon filters eliminate outside contamination.

## **Stand Alone or Integrated:**

The Standalone SuperVap® is ideal for manual, Semi-Automated, and Automated systems with no extract delivery capability. The Integrated SuperVap® is used with FMS, Inc. Fully automated Solid Phase Extraction Systems for PFAS analysis. The extract is automatically delivered to the centrifuge tube and concentrated or evaporated to a final volume without human intervention.

## **Industries**

Agricultural  
Clinical  
Environmental  
Food and Beverage  
Pharmaceutical Products  
Natural Products

# SuperVap® PFC Concentrator designed for PFAS Analysis

## Principals of Operation

SuperVap® Automated Concentration System is designed to streamline and automate the sample evaporation/concentration process for PFAS Analysis. It replaces techniques like KD, manual nitrogen blow-downs, and water baths. It automates existing manual evaporation and concentration processes in the laboratory. The design is intended to simplify, improve, and increase laboratory productivity. Automation of the Concentration evaporation process helps lower labor costs by minimizing the need for manual intervention. The SuperVap® brings efficiency, accuracy, and automation to the sample concentration process, reducing labor costs, saving time, and increasing overall productivity, making it a valuable addition to your laboratory workflow.

Compatible with existing FMS Sample Prep Systems

HEPA/Carbon Filter

HEPA/Carbon Filter to Eliminate Outside Contamination

Easy to Use Touch Screen Programming

Programmable Heat Ramp and Nitrogen Settings to Precisely Control the Concentration and Evaporation Process

A temperature log is saved for each run and may be downloaded to a PC via a USB port

Uses No Water, Dry heating assembly

Concentrates 12 to 24 Samples  
Sample Sizes: 15ml & 50ml

Compact Size

### Specifications

#### SuperVap®

Dimensions:

13"W x 13"D x 12"H"

Weight: 20 lbs.

Gas Requirements:

Nitrogen - 40 PSI minimum

Electrical Input:

110/220 Volts, 50/60 HZ

Controller: Touch Screen

Bath: Dry





# SuperVap® PFC Concentrator

## Automated Concentration and Evaporation

### EPA Method 533 Results

| Analyte        | Average Recoveries (%) | RSDs (%) |
|----------------|------------------------|----------|
| 11Cl-PF3OUdS   | 104.00                 | 3.50     |
| 9Cl-PF3ONS     | 86.90                  | 3.00     |
| ADONA          | 80.20                  | 2.90     |
| HFPO-DA (GenX) | 78.70                  | 3.70     |
| NFDHA          | 97.00                  | 5.00     |
| PFBA           | 101.20                 | 6.20     |
| PFBS           | 86.40                  | 3.70     |
| 8:2FTS         | 96.00                  | 4.00     |
| PFDA           | 89.50                  | 4.00     |
| PFDoA          | 80.60                  | 8.90     |
| PFEESA         | 98.00                  | 4.00     |
| PFHpS          | 100.20                 | 7.50     |
| PFHpA          | 83.00                  | 3.50     |
| 4:2FTS         | 98.00                  | 5.00     |
| PFHxS          | 92.00                  | 1.90     |
| PFHxA          | 101.20                 | 5.00     |
| PFMPA          | 98.00                  | 2.00     |
| PFMBA          | 101.00                 | 3.00     |
| PFNA           | 81.30                  | 3.00     |
| 6:2FTS         | 99.00                  | 11.00    |
| PFOS           | 86.20                  | 2.00     |
| PFOA           | 85.70                  | 2.90     |
| PFPeA          | 104.30                 | 6.30     |
| PFPeS          | 92.80                  | 5.50     |
| PFUnA          | 101.60                 | 9.20     |

### EPA Method 537.1 Results

| Compound      | Native Recoveries (%) |        |        |        |            |
|---------------|-----------------------|--------|--------|--------|------------|
|               | 2 ppt                 | 5 ppt  | 25 ppt | 50 ppt | EPA Window |
| PFBS          | 94.00                 | 93.00  | 98.00  | 99.00  | 70-130     |
| PFHxA         | 99.00                 | 104.00 | 101.00 | 109.00 | 70-130     |
| GenX          | 102.00                | 98.00  | 106.00 | 103.00 | NA         |
| PFHpA         | 99.00                 | 103.00 | 102.00 | 103.00 | 70-130     |
| PFHxS         | 95.00                 | 97.00  | 97.00  | 102.00 | 70-130     |
| ADONA         | 90.00                 | 97.00  | 99.00  | 104.00 | 70-130     |
| PFOA          | 116.00                | 109.00 | 105.00 | 103.00 | 70-130     |
| PFNA          | 95.00                 | 107.00 | 111.00 | 110.00 | 70-130     |
| PFOS          | 93.00                 | 96.00  | 95.00  | 101.00 | 70-130     |
| 9Cl-PF3ONS    | 88.00                 | 88.00  | 95.00  | 100.00 | 70-130     |
| PFDA          | 91.00                 | 99.00  | 105.00 | 111.00 | 70-130     |
| N-MeFOSAA     | 93.00                 | 97.00  | 92.00  | 92.00  | 70-130     |
| PFUdA         | 93.00                 | 101.00 | 104.00 | 108.00 | 70-130     |
| N-EtFOSAA     | 95.00                 | 110.00 | 98.00  | 98.00  | 70-130     |
| 11Cl-PF3OUdS  | 86.00                 | 88.00  | 86.00  | 91.00  | 70-130     |
| PFDoA         | 90.00                 | 92.00  | 99.00  | 101.00 | 70-130     |
| PFTeDA        | 86.00                 | 89.00  | 97.00  | 93.00  | 70-130     |
| PFTeDA (PFTA) | 84.00                 | 82.00  | 91.00  | 93.00  | 70-130     |

### EPA Method 1633 Results

| Compound     | Average Recoveries (%) | RSDs (%) | EPA Window (%) |
|--------------|------------------------|----------|----------------|
| PFBA         | 95.70                  | 1.00     | 70-135         |
| PFPeA        | 95.13                  | 2.00     | 70-135         |
| PFHxA        | 96.45                  | 3.00     | 70-135         |
| PFHpA        | 94.82                  | 1.00     | 70-135         |
| PFOA         | 95.78                  | 3.00     | 65-155         |
| PFNA         | 97.74                  | 2.00     | 70-140         |
| PFDA         | 97.69                  | 2.00     | 65-140         |
| PFUnA        | 96.42                  | 1.00     | 70-135         |
| PFDoA        | 96.65                  | 3.00     | 70-130         |
| PFTeDA       | 88.95                  | 4.00     | 60-145         |
| PFTeDA       | 98.78                  | 2.00     | 70-145         |
| PFBS         | 98.19                  | 3.00     | 70-140         |
| PFPeS        | 98.97                  | 1.00     | 70-135         |
| PFHxS        | 102.74                 | 2.00     | 70-135         |
| PFHpS        | 100.23                 | 6.00     | 70-140         |
| PFOS         | 100.65                 | 4.00     | 70-140         |
| PFNS         | 97.40                  | 3.00     | 70-135         |
| PFDS         | 88.92                  | 4.00     | 70-135         |
| PFDoS        | 84.00                  | 8.00     | 45-135         |
| 4:2FTS       | 95.46                  | 1.00     | 70-135         |
| 6:2FTS       | 98.17                  | 4.00     | 70-135         |
| 8:2FTS       | 92.35                  | 4.00     | 70-140         |
| PFOSA        | 98.99                  | 2.00     | 70-135         |
| NMeFOSA      | 90.47                  | 9.00     | 70-135         |
| NEtFOSA      | 92.16                  | 2.00     | 70-130         |
| NMeFOSAA     | 97.71                  | 4.00     | 65-140         |
| NEtFOSAA     | 97.11                  | 5.00     | 70-135         |
| NMeFOSE      | 92.73                  | 6.00     | 70-135         |
| NEtFOSE      | 91.80                  | 6.00     | 70-130         |
| HFPO-DA      | 99.24                  | 5.00     | 70-135         |
| ADONA        | 97.21                  | 5.00     | 70-135         |
| PFMPA        | 95.82                  | 1.00     | 60-140         |
| PFMBA        | 95.11                  | 2.00     | 65-145         |
| NFDHA        | 94.38                  | 6.00     | 65-145         |
| 9Cl-PF3ONS   | 95.92                  | 4.00     | 70-145         |
| 11Cl-PF3OUdS | 92.00                  | 3.00     | 50-150         |
| PFEESA       | 94.50                  | 3.00     | 70-135         |
| 3:3 FTCA     | 101.32                 | 2.00     | 70-130         |
| 5:3 FTCA     | 96.70                  | 3.00     | 70-130         |
| 7:3 FTCA     | 95.44                  | 3.00     | 55-130         |

# Total Solution Sample Prep for the Analysis of PFAS



**EconoTrace® PFC** for Drinking Water,  
Expandable from 2 to 8 Samples,  
Direct to Concentrator Delivery

**TurboTrace® PFC** for Drinking Water and  
Waste Water, Expandable for 5 to 30 Samples,  
Direct to Concentrator



## Applications

### For the analysis of


- Agricultural
- Clinical
- Environmental
- Food and Beverage
- Pharmaceutical Products
- Natural Products

## Automated Sample Preparation





# Semi-Automated Solid Phase Extraction for PFAS Analysis



The EZPFC-12<sup>®</sup> is a low-cost, high throughput semi-automated solid phase extraction solution for the analysis of PFAS.

## Systems

| Part Number        | Description   |
|--------------------|---|
| SVAP-PFC-12        | SuperVap <sup>®</sup> PFC Concentrator: 12 Position, 50ml Centrifuge Tubes                            |
| SVAP-PFC-24        | SuperVap <sup>®</sup> PFC Concentrator: 24 Position, 15ml Centrifuge Tubes                            |
| EZPFC-6            | 6 Position EZPFC System for Water & Waste Water Analysis  |
| EZPFC-12           | 12 Position EZPFC System for Water & Waste Water Analysis   |
| TT-SEQ-SPE-PFC     | TurboTrace <sup>®</sup> Parallel Sequential PFC Solid Phase Extraction System for 5 Samples           |
| TT-SEQ-SPE-PFC-EXP | TurboTrace <sup>®</sup> Parallel Sequential PFC Solid Phase Extraction Expansion Module for 5 Samples |
| ECO-PFC-SPE        | EconoTrace <sup>®</sup> PFC SPE Parallel System for 2 Samples   |
| ECO-PFC-SPE-EXP    | EconoTrace <sup>®</sup> PFC SPE Expansion Module for 2 Samples  |

## Accessories & Consumables

| Part Number  | Description  |
|--------------|--|
| HPCR-FIL-200 | Hepa/Carbon Filter set-complete                                      |
| SVAP-EXH-TUB | SuperVap <sup>®</sup> Polyethylene Exhaust tube, per foot            |
| FMS-TR-5012  | 12 Position Concentrator Tube Rack for 50ml Tube                     |
| VAC-PMP      | Vaccum Pump  |
| FMS_00773    | Fluid Line Kit EZSPE6 ESD 30", Sample Bottle to Cartridge, Pack of 6 |
| FMS_00635    | Fluid Line Kit PFC-12 30", Sample Bottle to Cartridge, Pack of 12    |
| SPE-ADP-SM   | SPE Small 1, 3, 6 ml Cartridge Adapter - Male to Male                |
| EZP-OWY-SCV  | One way stopcock valve   |



900 Technology Park Dr. Billerica, MA 01282  
www.fms-inc.com

Ph 617-393-2396  
Fax 617-393-0194  
sales@fms-inc.com